



ProMedic Nitrile Examination Gloves By Forte Global



Table of Contents

Product Specification

O2 Certificates

O3 SGS analysis

O4 Container Loading

Order Procedure & Forte Service Provided

06 Pictures



Product Specification

01.

M - ProMedic Nitrile Examination Gloves are made from 100% Synthetic Nitrile Latex which offers comfortability and excellent fit making it suitable for a wide range of applications. It provides excellent protection and is more durable than standard gloves. This is a suitable type of glove for many toxic tasks and works well with oil, grease, and even gasoline.

Quality Standards

- Conforms to ASTM D6319 and EN455 Standards
- Manufactured under QSR (GMP), ISO 9001:2015 and ISO 13485:2016 Quality Management System

Glove Sizes

- Extra-Small, Small, Medium, Large, Extra-Large
- Size of gloves shall be marked in the checkbox on the shipping carton with black ink



Product Specification

01.

Type

Powder-Free, Non-Sterile

Material

100% Synthetic Nitrile Latex
An alternative solution for those with a Latex Allergy

Color

Blue, White, Green, Pink, Light Purple, Black, Red

Packaging

100 pcs by weight per box, 10 boxes per carton



Product Specification

01.

Design & Features

- Protection from unwanted and dangerous substances
- Ambidextrous smooth material for hours of comfort
- Textured fingertips for a non-slip grip
- Beaded cuff for strength and easy donning

Storage

Store in a cool dry and location. Keep out of direct sunlight.

Shelf-Life

5 years from the date of manufacturing



ASTM D-6319

| Nitrile (Carboxylated Butadiene A Beaded Blue Ambidextrous | is a disposable of a disposable of and or finge on and examiner sion and Custome | device intend r to prevent | ed for medical purpos | | | | | | | |
|--|---|--|---|--|--|--|--|--|--|--|
| A patient examination glove that is worn on the examiner contamination between patie ASTM D-6319-00A Current Revis Nitrile (Carboxylated Butadiene ABeaded Blue Ambidextrous | is a disposable of a disposable of and or finge on and examiner sion and Custome | device intend r to prevent | ed for medical purpos | | | | | | | |
| that is worn on the examiner contamination between patie ASTM D-6319-00A Current Revis Nitrile (Carboxylated Butadiene A Beaded Blue Ambidextrous | 's hand or finge nt and examine sion and Custome | r to prevent | ed for medical purpos | | | | | | | |
| contamination between patie ASTM D-6319-00A Current Revis Nitrile (Carboxylated Butadiene A Beaded Blue Ambidextrous | nt and examine sion and Custome | | | | | | | | | |
| ASTM D-6319-00A Current Revis Nitrile (Carboxylated Butadiene A Beaded Blue Ambidextrous | sion and Custome | | | | | | | | | |
| Nitrile (Carboxylated Butadiene A Beaded Blue Ambidextrous | | Requirement | | | | | | | | |
| Beaded Blue Ambidextrous | Acrylonitrile) | | ASTM D-6319-00A Current Revision and Customer Requirement | | | | | | | |
| Blue Ambidextrous | | Nitrile (Carboxylated Butadiene Acrylonitrile) | | | | | | | | |
| Ambidextrous | | Beaded | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| Finger textured | | | | | | | | | | |
| | | | | | | | | | | |
| Size | Palm Width (n | nm) | Length (mm) | | | | | | | |
| X-Small | 75 ± 5mm | | 240min | | | | | | | |
| Small | 85 ± 5mm | | 240min | | | | | | | |
| Medium | 95 ± 5mm | | 240min | | | | | | | |
| Large | 106 ± 5mn | ı | 240min | | | | | | | |
| Extra Large | 115 ± 5mn | 1 | 240min | | | | | | | |
| Single wall | | | | | | | | | | |
| | Specificat | Specification | | | | | | | | |
| Cuff | | | 8 | | | | | | | |
| Palm | | | 0.06-0.09 | | | | | | | |
| Finger | 0.08-0.1 | 0.08-0.14 | | | | | | | | |
| Beading | 0.5min | 0.5min | | | | | | | | |
| | Un-aged | | Aged | | | | | | | |
| Tensile Strength MPA Min | 14 | | 14 | | | | | | | |
| | 500 | | 400 | | | | | | | |
| | | | | | | | | | | |
| Inspection Criteria | | Sa | mpling Plan/AQL | | | | | | | |
| Dimension | ' | | | | | | | | | |
| Force at Break | | 13 pcs (Median) | | | | | | | | |
| | | G1 – AQL 1.5 | | | | | | | | |
| | | G1 – AQL 2.5 | | | | | | | | |
| | G1 – | G1 – AQL 4.0 | | | | | | | | |
| Residual Powder | n = 5 | (average) | | | | | | | | |
| | | | | | | | | | | |
| | | n = 5 (AND) | | | | | | | | |
| | | | () | | | | | | | |
| 100 Gloves per Inner Dispense | r (by weight) 10 Di | spensers per C | uter Carton | | | | | | | |
| The dispensers and outer cartons shall have the following information/prints: | | | | | | | | | | |
| Product description and Brand including size and net quantity of gloves. | | | | | | | | | | |
| Country of Origin of Manufacturer and name and address of distributor. Symb | | | | | | | | | | |
| | | | • | | | | | | | |
| Per the customer's requiren | nent. | | | | | | | | | |
| _ | | nlicable | | | | | | | | |
| - Autoration and or expir | , auto miere ap | | | | | | | | | |
| | | | | | | | | | | |
| | | HUNGOUL | | | | | | | | |
| | X-Small Small Medium Large Extra Large Single wall Cuff Palm Finger Beading Tensile Strength MPA Min Elongation 2.0 mg/glc Inspection Criteria Dimension Force at Break Water Tight Test Visual Inspection – Major Visual Inspection – Minor Residual Powder Packaging Appearance Glove Physical Counting 100 Gloves per Inner Dispense The dispensers and outer ce Product description and Bra Country of Origin of Manu for single use and/or Single Per the customer's requiren | X-Small X-S | X-Small | | | | | | | |

EN 455

| Spec No. AQL12/02 | PRODUCT SPECIFICATION Revision: 3 Effective: 13/11/2020 | | | | | | | | |
|--------------------------------------|---|--|------------------|---------------------------------|--|--|--|--|--|
| Product Information: | I | | 2,1000 | | | | | | |
| Trade Name | Non-sterile Powder Free Nitrile E | xamination Glov | /es | | | | | | |
| Intended Use | A patient examination glove | is a disposable | e device intend | led for medical purpos | | | | | |
| | | that is worn on the examiner's hand or finger to prevent contamination between patient and examiner | | | | | | | |
| Product Conformance | EN 455 Current Revision and Cu | EN 455 Current Revision and Customer Requirement | | | | | | | |
| Material | Nitrile (Carboxylated Butadiene A | | | | | | | | |
| Cuff Finishing | Beaded | , | | | | | | | |
| Color | Blue | | | | | | | | |
| Туре | Ambidextrous | | | | | | | | |
| Surface finish | Finger textured | | | | | | | | |
| Product Specification: | | | | | | | | | |
| Dimension | Size | Palm Width | (mm) | Length (mm) | | | | | |
| | X-Small | 75 ± 5m | m | 240min | | | | | |
| | Small | 85 ± 5m | m | 240min | | | | | |
| | Medium | 95 ± 5m | m | 240min | | | | | |
| | Large | 106 ± 5r | nm | 240min | | | | | |
| | Extra Large | 115 ± 5r | nm | 240min | | | | | |
| Thickness | Single wall | Single wall thickness (mm) | | | | | | | |
| | | | tion | | | | | | |
| | Cuff | | 08 | | | | | | |
| | Palm | Palm | | | | | | | |
| | Finger | Finger | | | | | | | |
| | Beading | Beading | | | | | | | |
| Physical Properties | | Un-aged | | Aged | | | | | |
| | Force at break (N) | 6min | | 6min | | | | | |
| Residual Powder | 2.0 mg/gld | ove maximum | | | | | | | |
| Quality Specification | Inspection Criteria | | Sa | ampling Plan/AQL | | | | | |
| * AND – Accept No Defects | Dimension | | 13 pc | 13 pcs (Median) | | | | | |
| | Force at Break | | | 13 pcs (Median) G1 – AQL 1.5 | | | | | |
| | Water Tight Test | Water Tight Test | | | | | | | |
| | Visual Inspection – Major | | | | | | | | |
| | Visual Inspection – Minor | Visual Inspection – Minor | | | | | | | |
| | Residual Powder | Residual Powder | | | | | | | |
| | Packaging Appearance | Packaging Appearance | | | | | | | |
| | Glove Physical Counting | Glove Physical Counting | | | | | | | |
| Packaging Configuration and Labeling | | | | | | | | | |
| Packing | 100 Gloves per Inner Dispense | r (by weight) 10 | Dispensers per (| Outer Carton | | | | | |
| Labeling | The dispensers and outer cartons shall have the following information/prints: | | | | | | | | |
| | Product description and Brand including size and net quantity of gloves. | | | | | | | | |
| | Country of Origin of Manufacturer and name and address of distributor. Symbol | | | | | | | | |
| | for single use and/or Single | Use Only and | l Non-Sterile | | | | | | |
| Lot Numbering System | Per the customer's requiren Manufacturing and/or expir | | annlicable | | | | | | |
| Authorized by: | Code: | | | | | | | | |
| Tan Choong Kai QA Specialist | | | NBR35BL | | | | | | |
| Tan one of the art opening | | | HDIGGODE | | | | | | |
| | | | | | | | | | |
| | | | | | | | | | |



Certificates



















SGS analysis





JOB REF NO.: 2020-DATE RECEIVED: AUGUST 19, 2020 DATE REPORTED: SEPTEMBER 04, 2020 PAGE: 2 of 5

Test Report No. Company

TEST REPORT

Blue Nitrile Powder Free Examination Gloves Sample Description

Brand Name M-Promedic

Style

Lot No. 200818212

Colour Production Date August 2020

Malaysia Country of Origin

Size Quantity Tested 13 pieces

Test Conducted Dimensions

Test Method EN 455 Part 2:2015

Testing Period 19 August 2020 - 04 September 2020

Based on submitted samples, the following results obtained :-

| Size | M | M | M | M | M | M | M | M | M | M | M | M | M | Median |
|---------------------------|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|--------|
| Width Median: 95±10mm | 95 | 95 | 94 | 94 | 94 | 94 | 93 | 95 | 93 | 93 | 95 | 93 | 93 | 94 |
| Length Median: ≥ 240mm | 256 | 248 | 246 | 243 | 252 | 246 | 245 | 249 | 246 | 246 | 249 | 248 | 245 | 246 |

SIGNED FOR AND ON BEHALF OF SGS (MALAYSIA) SDN BHD

CHEE TUCK CHOON SECTION HEAD IKM No. M/3983/6401/12/14

content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 14 days only for perishable food samples, 30 days only for non-perishable food samples, and 90 days only for fats, oils and feed samples.

SGS (Malaysia) Sdn.Bhd. Lot 4, Persiaran Jubi Perak, Seksyen 22, 40300 Shah Alam, Selangor Darul Ehsan, Malaysia. (Company No. 10871-T) t+6(03) 7627 0080 f+6 (03) 7627 0082 www.sgs.com

Member of the SGS Group (SGS SA)



JOB REF NO.: 2020-0

DATE RECEIVED: AUGUST 19, 2020 DATE REPORTED: SEPTEMBER 04, 2020

PAGE: 5 of 5

Test Report No. Company

TEST REPORT

Sample Description Blue Nitrile Powder Free Examination Gloves

Brand Name M-Promedic

M Style

200818212 Lot No. Colour Blue

Production Date August 2020 Country of Origin Malaysia

M Size

Quantity Tested 8 pieces

Test Conducted Aqueous Extractable Protein Content

Test Method EN455 Part 3:2015

19 August 2020 - 04 September 2020 Testing Period

EXTRACTABLE PROTEIN CONTENT

Analysis had been carried out as per your request. We report the following results:

| SIZE | AQUEOUS EXTRACTABLE PROTEIN CONTENT, μg/g |
|------|---|
| M | 56 |

Note: Detection limit = 10 µg/g

SIGNED FOR AND ON BEHALF OF SGS (MALAYSIA) SDN BHD

CHEE TUCK CHOON SECTION HEAD

IKM No. M/3983/6401/12/14

This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com/en/Terms.and-Conditions/Terms.and-Condit confert or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 14 days only for perishable food samples, 30 days only for non-perishable food samples, and 90 days only for fats, oils and feed samples.

SGS (Malaysia) Sdn.Bhd. Lot 4, Persiaran Jubli Perak, Seksyen 22, 40300 Shah Alam, Selangor Darul Ehsan, Malaysia

(Company No. 10871-1) t+6(03) 7627 0080 f+6 (03) 7627 0082 www.sgs.com





JOB REF NO.: 2020-08-

DATE RECEIVED: AUGUST 19, 2020 DATE REPORTED: SEPTEMBER 04, 2020

Test Report No. Company

TEST REPORT

Sample Description Blue Nitrile Powder Free Examination Gloves

M-Promedic Brand Name Style

Lot No. 200818212 Colour Blue Production Date August 2020 Country of Origin Malaysia M

13 pieces Quantity Tested

Test Conducted Dimensions ASTM D6319-10 (Reapproved 2015) Test Method Testing Period 19 August 2020 - 04 September 2020

Based on submitted samples, the following results obtained :

| Dascu on submitted sam | pies, the | OHOWI | ng resu | ns oou | med | | | | | | | | |
|------------------------------------|-----------|-------|---------|--------|------|------|------|------|------|------|------|------|------|
| Size | M | M | M | M | M | M | M | M | M | M | M | M | M |
| Width 95 ± 10mm | 93 | 93 | 94 | 95 | 94 | 93 | 94 | 94 | 93 | 93 | 93 | 93 | 93 |
| Length Min. 230mm | 246 | 247 | 253 | 249 | 245 | 245 | 250 | 251 | 250 | 244 | 251 | 247 | 244 |
| Thickness at palm Min. 0.05mm | 0.07 | 0.07 | 0.07 | 0.06 | 0.07 | 0.08 | 0.07 | 0.06 | 0.07 | 0.07 | 0.07 | 0.07 | 0.07 |
| Thickness at finger Min. 0.05mm | 0.08 | 0.09 | 0.08 | 0.08 | 0.09 | 0.10 | 0.08 | 0.08 | 0.09 | 0.08 | 0.08 | 0.09 | 0.08 |

Acceptable Quality Limit (AQL): 4.0

Accept: 1

Found: 0

SIGNED FOR AND ON BEHALF OF SGS (MALAYSIA) SDN BHD

CHEE TUCK CHOON SECTION HEAD

IKM No. M/3983/6401/12/14

content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 14 days only for perishable food samples, 30 days only for non-perishable food samples, and 90 days only for fats, oils and feed samples.

SGS (Malaysia) Sdn.Bhd. Lot 4, Persiaran Jubii Perak, Seksyen 22, 40300 Shah Alam, Selangor Darul Ehsan, Malaysia. (Company No. 10871-7) 1+6(03) 7627 0080 1+6 (03) 7627 0082 www.sgs.com

Member of the SGS Group (SGS SA)



JOB REF NO.: 2020-

DATE RECEIVED: AUGUST 19, 2020 DATE REPORTED: SEPTEMBER 04, 2020

PAGE: 3 of 5

Test Report No. Company

TEST REPORT

Sample Description Blue Nitrile Powder Free Examination Gloves

Brand Name M-Promedic

Style

Lot No. 200818212 Colour Production Date August 2020 Malaysia

Country of Origin

Quantity Tested 13 pieces per each

Test Conducted Force at Break During Shelf Life and After Challenge

Test Method EN 455 Part 2:2015 70 ± 2 Deg C for 168 hrs

Testing Period 19 August 2020 - 04 September 2020

| | | Force at | Break, N |
|-------------|------------|--------------|-------------|
| SIZE | SAMPLE NO. | BEFORE AGING | AFTER AGING |
| M | 1 | 6.7 | 6.3 |
| | 2 | 6.3 | 6.5 |
| | 3 | 6.9 | 6.6 |
| | 4 | 6.8 | 7.6 |
| | 5 | 8.2 | 6.5 |
| | 6 | 6.2 | 6.8 |
| | 7 | 6.1 | 7.4 |
| | 8 | 7.3 | 6.1 |
| | 9 | 6.0 | 7.6 |
| | 10 | 7.0 | 8.7 |
| | 11 | 6.2 | 6.5 |
| | 12 | 7.6 | 6.3 |
| | 13 | 6.0 | 8.1 |
| Median | | 6.7 | 6.6 |
| Requirement | | >60 | >60 |

SIGNED FOR AND ON BEHALF OF SGS (MALAYSIA) SDN BHD

CHEE TUCK CHOON SECTION HEAD

IKM No. M/3983/6401/12/14

This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.aga.com/wnTerms.and-Conditions.tom and Conditions for Bestoreic Documents as <a href="http://www.aga.com/wnTerms.and-Conditions/Terms. obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 14 days only for perishable food samples, 30 days only for non-perishable food samples, and 90 days only for fats, oils and feed samples.

SGS (Malaysia) Sdn Bhd. Lot 4, Persiaran Jubli Perak, Seksyen 22, 40300 Shah Alam, Selangor Darul Ehsan, Malaysia. (Company No. 10871-T) t+6(03) 7627 0080 f+6 (03) 7627 0082 www.sqs.com





JOB REF NO.: 2020-DATE RECEIVED: AUGUST 19, 2020 DATE REPORTED: SEPTEMBER 04, 2020

Test Report No. Company

TEST REPORT

Sample Description Blue Nitrile Powder Free Examination Gloves

Brand Name M-Promedic

Style

Lot No. 200818212

Colour Blue Production Date August 2020

Country of Origin Malaysia

Size

Quantity Tested 200 pieces

Test Conducted Freedom from holes

Test Method EN455 Part 1:2000

19 August 2020 - 04 September 2020 Testing Period

Based on submitted samples, the following results obtained :-

Acceptable Quality Limit (AQL): 1.5 Accept: 7 Found: 2

Result : Within AOL

SIGNED FOR AND ON BEHALF OF SGS (MALAYSIA) SDN BHD

CHEE TUCK CHOON SECTION HEAD

IKM No. M/3983/6401/12/14

SGS (Malaysia) Sdn Bhd. | Lot 4, Persiaran Jubli Perak, Seksyen 22, 40300 Shah Alam, Selangor Darul Ehsan, Malaysia. (Company No. 10871-T) t +6(03) 7627 0080 f +6 (03) 7627 0082 www.sqs.com

Member of the SGS Group (SGS SA)



JOB REF NO.: 2020-08-

DATE RECEIVED: AUGUST 19, 2020 DATE REPORTED: SEPTEMBER 04, 2020

PAGE: 1 of 5

Test Report No. Company

TEST REPORT

Sample Description Blue Nitrile Powder Free Examination Gloves

Brand Name M-Promedic

Lot No. 200818212 Colour Blue

Production Date August 2020 Country of Origin Malaysia Size

Quantity Tested 200 pieces Test Conducted Freedom from holes

Test Method ASTM D6319-10 (Reapproved 2015)

Testing Period 19 August 2020 - 04 September 2020

Based on submitted samples, the following results obtained :-

Acceptable Quality Limit (AQL): 2.5 Accept: 10 Found: 2

: Within AQL

SIGNED FOR AND ON BEHALF OF SGS (MALAYSIA) SDN BHD

CHEE TUCK CHOON SECTION HEAD IKM No. M/3983/6401/12/14

This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <a href="https://document.gov/rene-printed-conditions-page-and, for electronic format documents, subject to Tames and Conditions for Electronic Documents at <a href="https://documents.gov/rene-printed-conditions-format-conditions-fo Isability, indemnification and particulities issues defined therein. Any holder of this document is activated that information contained between reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole regionability is to its Client and this document does not exceed parties to a transaction from exercizing all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgety or faitafication of the content or appearance of this document is unlessful and effectives may be presecuted to the fullest estent of the law. Unless otherwise stated the results shown in this test report refer only to the amplicity listed and such sample(s) are retained for 1 days only for perinhable food samples, 30 days only for non-perishable food samples, 30 days only for non-perishable food samples. 30 days only for non-perishable food samples, 30 days only for non-perishable food samples, 30 days only for non-perishable food samples. 30 days only for non-perishable food samples, 30 days only for non-perishable food samples. 30 days only for non-perishable food samples. 30 days only for non-perishable food samples, 30 days only for non-perishable food samples. 30 days only for non-perishab

SGS (Malaysia) Sdn.Bhd.
(Company No. 10871-T)

Lot 4, Persiaran Jubii Perak, Seksyen 22, 40300 Shah Alam, Selangor Darul Ehsan, Malaysia.

1+6(03) 7627 0080 #+6 (03) 7627 0082 www.sgs.com





JOB REF NO.: 2020-08-DATE RECEIVED: AUGUST 19, 2020 DATE REPORTED: SEPTEMBER 04, 2020

PAGE: 3 of 5

Test Report No. Company

TEST REPORT

Blue Nitrile Powder Free Examination Gloves Sample Description

Brand Name M-Promedic

Style

Lot No. 200818212 Colour Blue August 2020 Production Date Country of Origin Malaysia

M

Sample Quantity 13 pieces per each

Test Conducted Tensile Strength & Elongation (Before Ageing & After Ageing)

Test Method ASTM D6319-10 (Reapproved 2015)

Ageing 70 ± 2 Deg C for 168 hrs 19 August 2020 - 04 September 2020 Testing Period

| Size | Sample | BEFORE | AGEING | AFTER | AGEING |
|--------------|--------|----------|------------|----------|------------|
| | No. | Tensile | Ultimate | Tensile | Ultimate |
| | | Strength | Elongation | Strength | Elongation |
| | | (MPa) | (%) | (MPa) | (%) |
| M | 1 | 24.2 | 520 | 25.3 | 480 |
| | 2 | 23.8 | 520 | 24.7 | 480 |
| | 3 | 21.8 | 500 | 22.6 | 460 |
| | 4 | 30.2 | 560 | 25.4 | 480 |
| | 5 | 27.7 | 540 | 27.2 | 480 |
| | 6 | 24.9 | 520 | 29.6 | 500 |
| | 7 | 27.9 | 540 | 30.4 | 480 |
| | 8 | 28.3 | 540 | 32.1 | 500 |
| | 9 | 26.6 | 540 | 32.3 | 520 |
| | 10 | 31.3 | 560 | 26.5 | 460 |
| | 11 | 31.0 | 560 | 27.5 | 480 |
| | 12 | 24.4 | 500 | 26.6 | 460 |
| | 13 | 24.5 | 500 | 26.6 | 460 |
| Requirements | | 14 min | 500 min | 14 min | 400 min |

SIGNED FOR AND ON BEHALF OF SGS (MALAYSIA) SDN BHD

CHEE TUCK CHOON SECTION HEAD IKM No. M/3983/6401/12/14

This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessable at http://www.aps.com/em/Terms.and-Conditions.aps.and, for electronic format documents, subject to Terms and Conditions for Sections Documents as all the format and the section of the section of the section of the section and practication and practication is above to the third and the section of th

SGS (Malaysia) Sdn. Bhd. Lot 4, Persiaran Jubli Perak, Seksyen 22, 40300 Shah Alam, Selangor Darul Ehsan, Malaysia. (Company No. 10871-T) t +6(03) 7627 0080 f +6 (03) 7627 0082 www.sgs.com

Member of the SGS Group (SGS SA)



JOB REF NO.: 2020-

DATE RECEIVED: AUGUST 19, 2020 DATE REPORTED: SEPTEMBER 04, 2020

PAGE: 4 of 5

Test Report No. Company

TEST REPORT

Blue Nitrile Powder Free Examination Gloves Sample Description

Brand Name M-Promedic

Style Lot No. 200818212

Colour Production Date August 2020 Country of Origin Malaysia

Size Quantity Tested 5 pieces

Test Conducted Powder Content Test Method EN455 Part 3:2015

19 August 2020 - 04 September 2020 Testing Period

On testing the samples, the following results were obtained:-

SIZE Average Powder Mass per Glove 0.38 mg M

SIGNED FOR AND ON BEHALF OF SGS (MALAYSIA) SDN BHD

CHEE TUCK CHOON SECTION HEAD

IKM No. M/3983/6401/12/14

This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <a href="http://www.aga.com/wnTerms-and-Conditions-Imms-and-Condit obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 14 days only for perishable food samples, 30 days only for non-perishable food samples, and 90 days only for fats, oils and feed samples.

SGS (Malaysia) Sdn Bhd.
(Company No. 10871-7)

Lot 4, Persiaran Jubii Perak, Seksyen 22, 40300 Shah Alam, Selangor Darul Ehsan, Malaysia.

1, 46(03) 7627 0080 1+6 (03) 7627 0080 www.sgs.com




Container Loading



Dry Freight Container





20ft



Container 20 ft Standard:

11,880 Boxes (1,188 Cartons)

Container 40 ft Standard High Cube:

28,880 Boxes (2,880 Cartons)

Carton Dimensions:

390mm x 250mm x 250mm





Order Procedure & Forte Service Provided



Order Procedure & Forte Service Provided

05.

- 1. Buyer Issues Purchase Order (PO) to Forte
- 2. Forte Issues Proforma Invoice (PI) to Buyer
- Buyer sends 50% deposit of the PO to Forte within 48 hours. Once the
 deposit is received, the PO is officially locked into production and
 Forte will advise Buyer the estimated date the production will be
 ready by
- 4. Buyer arranges their own shipping. However, Forte can assist using their own Shipping Agent and bill them directly. The invoice will be for Ocean Freight, Misc. charges, and Export Rubber Tax (0.2% of the PO)
- 5. Buyer advises Forte when the container will arrive so the necessary arrangements can be made to prepare for the loading day.
- 6. Once stock is ready, Forte advises Buyer



Order Procedure & Forte Service Provided

05.

- 7. Buyer arranges their own 3rd party inspection (if needed). If Buyer proceeds with their own inspection, Forte will coordinate with the Buyer and the Inspection company to set the date
- 8. Forte will set a date before the official 3rd inspection to perform their own pre-inspection with Buyer on a live zoom call
- 9. Official 3rd inspection is done with Forte and the Inspection company. Buyer will have a live zoom call available for them to witness the entire inspection
- 10. Buyer sends remaining balance (50%) of the PO to Forte after inspection
- 11. On loading day, Forte will have a live zoom call available for Buyer to witness the entire loading and sealing of the container



Order Procedure & Forte Service Provided

05.

- 12. Depending on when the container is picked picked up, whether day or another date, Forte will have a live zoom call for Buyer to witness the container being picked up by the carrier and leave the facility to head to the port
- 13. Deal is then concluded. Forte will continue to assist in any way possible after this point and will also track the container as a courtesy to the Buyer all the way to the final destination

Note:

In any case where Buyer is not able to be present for the live zoom call, a video recording will be made available for Buyer to view at their own convenience



Pictures













ProMedic

Nitrile Examination Gloves

















































Thank you. Questions/Comments?

hermanliu@forteglb.com +60 12 775 4248 forteglb.com

